




Item/Color	Size	SWL
AM-TSL-L-F8	46" x 86"	1000 Lb/ 454 Kg
AM-TSL-W-F8	50" x 86"	1000 Lb/ 454 Kg
AM-TSL-TS-L-F8	46" x 86"	1000 Lb/ 454 Kg
AM-TSL-TS-W-F8	50" x 86"	1000 Lb/ 454 Kg


### SLING PLACEMENT ON BED

- PLACE BARRIER SHEET OR CHUCK BETWEEN PATIENT AND SLING - DO NOT LAY PATIENT DIRECTLY ON MESH.
- FOOT LABEL SHOULD BE AT FOOT END.
- PLACE PILLOW BETWEEN KNEES WHEN LOG ROLLING PATIENT FOR HIP PRECAUTIONS.




### LOG ROLLING

- TO LOGROLL : CHOOSE TWO STRAPS, ONE CLOSE TO PATIENT SHOULDER AND ONE CLOSE TO PATIENT HIP.
- USE SHORTEST ATTACHMENT POINT ON FIGURE 8 STRAP.



### SUPINE TRANSFER

- ATTACH **STRAPS 1-2** ON BOTH SIDES TO THE 2 HOOKS NEAREST TO PATIENT'S HEAD.
- ATTACH **STRAPS 3, 4 & 5** ON EITHER SIDE TO THE 2 HOOKS NEAREST TO FEET.



### RECLINED TRANSFER

- ATTACH **STRAPS 1-2** ON BOTH SIDES TO THE 2 HOOKS NEAREST TO PATIENT'S HEAD.
- ATTACH **STRAPS 3, & 4** ON EITHER SIDE TO THE 2 HOOKS NEAREST TO FEET.
- ATTACH **STRAP 5** IF PATIENT TOLERATES LEG HYPEREXTENSION - OTHERWISE LEAVE STRAP UNATTACHED.