

## A-BedLadderSW





SWL 600 LB / 272 KG

Cleaning:  
Wipe with facility  
approved disinfectant  
between patients  
**Do not launder**



## SUPINE TO SITTING

Attach Bed Ladder to the patient's bed frame:

- 1  Push the end straps through the opening of the bed frame. Pull straps up and over the outside of the openings.
- 2  Pull pronged end of clip through one strap and repeat with the other strap.
- 3  Ensure that the pronged end of the clip is going through BOTH end straps.
- 4  Connect the pronged end of clip with the open end of clip so that clip attachment is clasped shut.

Have patient use Bed Ladder to pull themselves up into a sitting position:

NOTE: Bed rails are down in pictures in order to show maneuvers. Please adhere to your hospital policy regarding bed rail safety.



A-BedLadderSW

SWL 600 LB / 272 KG

LIFTING LEGS BACK ONTO BED

1


Have patient place one foot into both loops of Bed Ladder.

2


Bed Ladder can be used to raise foot and swing it up onto the bed.

3


Have patient repeat the process with the other leg.



Bed Ladder can also be used to help patient get into a vehicle.

SITTING UP ON EXAM TABLE



Place the loops of the Bed Ladder around the patient's feet



Have patient use Bed Ladder to pull up to an upright sitting position

Disclaimer

1. Discontinue use if loose stitching or tears are visible  
2. DO NOT leave patient unattended with a Bed Ladder  
3. When using Bed Ladder, rails should be in the up position