

A-SitOnSling



ITEM	Seat Width	Recommended Weight Range	Max SWL
A-SitOnSling-M	30" / 76cm	75-145lbs / 34-66kg	1000 lbs / 454kg
A-SitOnSling-L	32" / 81cm	145-220lbs / 66-100kg	
A-SitOnSling-X	35" / 89cm	220-440lbs / 100-200kg	
A-SitOnSling-XX	39" / 99cm	440-800lbs / 200-363kg	

*Trim color indicates size

Placing SitOnSling under patient



Place sling on the bed with label side down.



Make sure straps on the side to be rolled are folded into the sling so that they do not touch the patient.
Roll up sling approximately one quarter of the way.



Roll patient to one side and place SitOnSling behind the patient between shoulders and tail bone. Make sure the label remains face down on the bed.



Roll patient over the sling and over to their other side.



Safely unroll sling under patient.



Roll patient onto their back.

Using A-SitOnSling to lift patient

With patient on their back, raise head of bed to a comfortable angle. Attach loops as follows:



4pt hanger bar: Attach each of the four straps (both shoulder and both leg) to each one of the hooks on the hanger bar.



2pt hanger bar: Attach the left side shoulder and leg straps to the left side hook. Attach the right side shoulder and leg straps to the right side hook.



Note: the lower the loop you use on the strap, the closer the patient will be to the hanger bar.

Typically, if you use a lower loop on the shoulder straps than the leg straps, the patient's head will be higher up.

Cleaning Instructions

- Wash Temp not to exceed 185F / 85C
- Drying - Medium Heat, not to exceed 165F / 74C
- Oxygen Based Bleach Recommended
- Max Load 850lbs Cycle