



Your Sling Source for Safe Patient Mobilization

New Sling Labels- Washable/Wipeable

October 2018

In order to be in compliance with the new guidelines set forth by the FDA regarding UDI barcodes, we have begun to implement our new label design. Initially, you will begin to see these new labels on our washable seated slings and washable TSL slings. Eventually, all slings, including custom slings, will have these new labels with the only exception being our SPU/ECO line.

Alpha Modalities
1111 S. 344th St., Ste 205
Federal Way, WA 98003
1-800-273-5749

A-SEAT Series
Attaches to a 2, 4 or 6 pt Loop Compatible Hanger Bar
Max Load 1000lbs / 454kg

! This sling must be inspected before each use
Do Not Use if Damaged - Do Not Repair

Made in USA

LARGE
Recommended Weight Range:
145-220 lbs / 66-100 kg

For more information on use, laundry or maintenance of this sling, go to www.alphamodalities.com

85C [Washing symbols]

	Date of First Use MM/YR		
--	----------------------------	--	--

Universal symbols for cleaning instruction

Indicates type of attachment (loop or clip) and Maximum Safe Working load

Indicates specific size of sling. Box color indicates trim color. Recommended weight for this size is also listed.

This box will have the UDI barcode inserted specific to the sling type, fabric and size.

