


Your Sling Source for Safe Patient Mobilization

1111 South 344th Street, Suite 205, Federal Way, WA 98003  
800.273.5749 • info@alphamodalities.com


### Pannus Single Patient Use Sling

Description	Part Number	Max Load	Order Quantity
Pannus Sling	A-PannusSPU	454kg / 1000lb	Sold as individual item

COMPATIBLE WITH LOOP SYSTEMS  
OR CAN BE USED MANUALLY



INSPECTION SLITS  
ON SLING FOR CLEANING  
AND/OR EXAM



**BREATHABLE MATERIAL**  
SOFT TO THE TOUCH  
LATEX FREE



**DO NOT LAUNDR**  
SPOT CLEAN WITH  
HOT WATER IF SOILED



**THROW SLING AWAY**  
AFTER PATIENT USE

