


Your Sling Source for Safe Patient Mobilization

1111 South 344th Street, Suite 205, Federal Way, WA 98003
800.273.5749 • info@alphamodalities.com

TSL Single Patient Use Repositioning Sling

Description	Part Number	Size	Max Load	Order Quantity
Repositioning Sling	A-TSLSPU	45"W x86"L	454kg / 1000lbs	Sold by the case: 10 slings/case
	A-TSLSPUX	55"W x86"L	500kg / 1100lbs	


MATTRESS STRAPS
KEEP SLING FROM SLIDING DOWN THE BED




2 LOOPS PER STRAP
PROVIDE FOR BETTER POSITIONING



THROW SLING AWAY
AFTER PATIENT USE



FITS LENGTH OF THE BED
86" LENGTH



BREATHABLE MATERIAL
SOFT TO THE TOUCH
LATEX FREE



PRESSURE MAPPING TESTING
SLING CAN REMAIN
UNDERNEATH PATIENT

